


System configuration



“SmartSwitch” Technology

Fujifilm developed a new technology “SmartSwitch” which allows automatic X-ray detection. With “SmartSwitch,” FDR D-EVO II no longer requires connection between the X-ray generator and DR power supply unit to automatically detect X-rays and start image creation.



Specification

| Model name | Flat Panel Detector (DR-ID 1201SE) for FDR D-EVO II System (DR-ID1200) | Flat Panel Detector (DR-ID 1202SE) for FDR D-EVO II System (DR-ID1200) |
|-------------------------|---|---|
| Type | Cassette size detector with ISS (Irradiation Side Sampling system) | Cassette size detector with ISS (Irradiation Side Sampling system) |
| Scintillator | GOS (Gadolinium oxysulfide) | GOS (Gadolinium oxysulfide) |
| Detector external size | 460 × 384 × 15 mm (Approx.) [18" × 15" × 0.6"] | 460 × 460 × 15 mm (Approx.) [18" × 18" × 0.6"] |
| Weight | Approx. 2.6kg [5.7lbs.] (including battery) | Approx. 3.2kg [7.1 lbs.] (including battery) |
| Pixel pitch | 0.15 mm | 0.15 mm |
| Pixels | 2836 × 2336 pixels | 2836 × 2832 pixels |
| Wireless standard | IEEE 802.11n (2.4GHz, W52/W53/W56/W58) | IEEE 802.11n (2.4GHz, W52/W53/W56/W58) |
| Image preview | Less than 2sec | Less than 2sec |
| Cycle time | Less than 9 sec (wired) Less than 10 sec (SmartSwitch) | Less than 9 sec (wired) Less than 10 sec (SmartSwitch) |
| Battery recharging time | Approx. 3 hours (with battery charger) Approx. 4 hours (with Docking Stand) | Approx. 3 hours (with battery charger) Approx. 4 hours (with Docking Stand) |
| Battery performance | Standby: Approx. 4 hours Sleep mode: Approx. 7h 30min Extra sleep mode: Approx. 18h 30min | Standby: Approx. 4 hours Sleep mode: Approx. 7h 30min Extra sleep mode: Approx. 18h 30min |

Optional parts



Battery charger Battery Fujifilm AP

External appearance and specifications are subject to change without notice.
 All brand names or trademarks are the property of their respective owners.
 All products require the regulatory approval of the importing country.
 For details on their availability, contact our local representative.
 Please contact FUJIFILM's authorized distributor for FDR D-EVO II X-ray system.




Be Smart.



NEW

FDR D-EVO II
G35 | G43



High definition, made smarter.



G35 [14"×17" model]

G43 [17"×17" model]


Improved image quality through a noise reduction circuit, and a variety of functions to support imaging. Introducing FDR D-EVO II, now even smarter.

- Designed to be lightweight; only 2.6 kg* with replaceable battery *14"×17" model
- Loaded with internal memory that allows detector-only image storage
- Antibacterial, waterproof, and load resistant performance for peace of mind during use
- LED indicators on detector edge confirm center location and distinguish multiple detectors in department
- The rounded form of the detector edges makes handling and patient positioning easy

Fujifilm's exclusive technology for achieving high resolution and low dosing


1. ISS system reading technology promotes high sensitivity

Equipped with an indirect conversion system flat detector display using ISS, which bonds optical sensors (TFT) to the X-ray irradiation side unlike traditional flat detector displays. This greatly suppresses scattering and attenuation of X-ray signals, creating sharp images with low doses of X-rays.



2. Blending large and small phosphor particles at an optimal ratio



Increased X-ray absorption through our uniquely designed scintillator, which applies photographic film technology to optimize the blending ratio of phosphor particles of different sizes.



Blending optimally-sized phosphor particles in gaps

3. Fujifilm noise reduction circuit improves detector sensitivity in high absorption regions

The uniquely developed noise reduction circuit reduces noise in the image. It achieves 1.7 times the DQE of existing systems with a 0.03 mR dose. In particular, granularity of low-concentration regions such as the heart and mediastinum is dramatically improved.



With additional major increases in sensitivity in low-concentration regions (heart, mediastinum)

4. Image processing technology to optimize imaging results

FDR D-EVO II utilizes the latest Fujifilm digital image processing technologies including Dynamic Visualization, which optimizes image display based on monitor characteristics and FNC noise suppression processing that improves image quality, automatically extracting and separating noise components in the image.

