



VILLA
SISTEMI
MEDICALI

MEDICAL
IMAGING


Apollo open

Open Tabletop



The new "open" tabletop for Apollo extends the application range to interventional procedures and provides total accessibility for improved patient positioning.

La nuova barella "open" estende l'ambito applicativo alle procedure interventistiche e offre la massima accessibilità per il posizionamento dei pazienti.



- Unrestricted patient access from all sides
- Patient weight capacity up to 230kg without movement limitations
- Fast, comfortable and safe patient access
- Interventional applications capability
- Available for Apollo and Apollo DRF
- Carbon fiber tabletop

- *Accesso totale al paziente da ogni lato*
- *Portata massima paziente fino a 230kg, senza limitazioni di movimento*
- *Accesso al paziente rapido, confortevole e sicuro*
- *Possibilità di procedure interventistiche*
- *Disponibile per Apollo e Apollo DRF*
- *Piano porta-paziente in fibra di carbonio*

Villa Sistemi Medicali: long-standing experience at the service of our customers

Villa Sistemi Medicali Spa
20090 Buccinasco (Mi) - Via delle Azalee, 3
Tel. +39 02 48859.1 - Fax +39 02 4881.844
www.villasm.com - vsminfo@villasm.com



Company with Quality System certified by:

