

VetPro[®] 1000

Mobile Dental
Delivery System



The VetPro® 1000 Mobile Dental Delivery System

Oil-Free, Worry-Free, Hassle-Free

A reliable veterinary dental delivery solution that helps you provide better patient care.



Oil-free compressor

Our exclusive oil-free compressor technology allows you to provide your patients with better, more efficient care without the mess or hassles of oil-lubricated units. The VetPro® 1000 dental delivery system requires no oil changes and very little maintenance for hassle-free operation. Oil-free means there is no oil to contaminate your handpieces, providing a clean, dry, aseptic work environment for you, your staff and your patients.

Designed for veterinary dentistry

The VetPro® 1000 has precise controls designed specifically for veterinary dentistry. You can switch from handpiece to handpiece without touching a button – the unit will automatically adjust as the handpiece is lifted from its cradle. The oil-free compressor means there are no oily parts for animal hair to stick to, and with a sound level of 40 dB, it's one of the quietest compressors on the market.

Durable construction

Its durable construction withstands the rigors often found in clinical settings. The VetPro® 1000 is a completely enclosed system to protect internal mechanics from contaminants and animal hair, and is made from sturdy ABS plastic to withstand the strains caused by frequent usage.



Find us on
Facebook

facebook.com/MidmarkAnimalHealth

midmarkanimalhealth.com

Midmark is an ISO 13485 and ISO 9001 Certified Company. • 1-800-MIDMARK • 1-937-526-3662 • midmark.com
© 2013 Midmark Corporation • Products subject to improvement changes without notice • Litho in U.S.A. • 007-1025-00 Rev. A1 (4/13)

Supported by Midmark

The VetPro® 1000 gives you unsurpassed reliability, backed by our trained service technicians and independent service network. You will have peace of mind knowing that over 90 years of manufacturing experience goes into each product. VetPro® 1000—the oil-free, worry-free, hassle-free solution for your veterinary dental needs.

Specifications

Overall depth: 18.5" (46.99 cm)

Overall width: 23" (58.42 cm)

Overall height: 36.5" (92.71 cm)

Casters: 3" (7.62 cm)
medical grade

Compressor

Oil-free 1/3 hp motor

Working pressure: 70-100 psi

Flow rate: 2.5 CFM at 40 psi

Sound level: 40 dB

Air reservoir: 1 gallon

Manual air tank purge

Water Module

Two 1-liter reservoirs (toggled)

Pressurized system: 30-40 psi

Aseptic, non-retracting
HPC coolant

Manual handpiece flush system

External quick connects for
water and air

Electrical

Voltage: 115 VAC, 6 Amps

Frequency: 50/60 Hz

UL 60601-1:2003

CAN/CSA-C22.2 No. 601-1-M90

Powercord: 10' hospital grade

Accessories

VetPro® Lube-Free High-Speed and
Low-Speed Handpieces included

Integrated Acteon Satelec®
SP Newtron® Piezo Scaler,
with or without LED Light

Handpiece connection type:
Midwest 4-hole style

Optional fiber optics

Bur starter kit included

10 Disposable prophylaxis angles
included

Optional antimicrobial fabric cover

Model Numbers

8000-001 VetPro® 1000
with integrated scaler

8000-003 VetPro® 1000
with integrated scaler
and fiber optics

8000-007 VetPro® 1000
with LED scaler and dual
fiber optic high-speed
handpieces

8000-008 VetPro® 1000
with LED scaler and
fiber optics



Because we care.