



## **DP-10Vet**

Digital Ultrasonic Diagnostic Imaging System

**mindray**  
healthcare within reach

# DP-10Vet

## Digital Ultrasonic Diagnostic Imaging System

### Ergonomic design with fashionable outlook

- 12.1" LED monitor with 30 degree tilting function
- New style backlit, silicon keyboard providing fantastic operations
- Magnet design for easy fold-up
- Slim system for extreme portability



### Excellent image quality and wide veterinary applications

- Full digital technology
- IP (Image Processing) for fast image optimization
- 8-TGC for accurate image control
- A broad range of transducers covering abdomen, thoracic, cardiac, tendinous, small parts, ophthalmic, reproduction system exams for animals such as dog, cat, equine, bovine and ovine, etc.

### Saving the scanning time by facilitating the workflow

- Fast boot-up within 20s
- User-defined keys for personalized operations
- One key images or cines saved to local disk or USB



DISTRIBUTOR:



MINDRAY is a trademark of Shenzhen Mindray Bio-Medical Electronics Co., Ltd. Specifications subject to changes without prior notice.  
© 2010-2011 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. All rights reserved.  
P/N:DP-10vet-ENG-210285-20110921

## mindray

Mindray is listed on the NYSE under the symbol "MR"  
Mindray Building, Keji 12th Road South, High-tech Industrial Park,  
Nanshan, Shenzhen 518057, P.R. China  
Tel: +86 755 8188 8998 Fax: +86 755 26582680  
E-mail: intl-market@mindray.com Website: www.mindray.com